

Certificate of Analysis for NR-51666

Yersinia pestis, Strain PY-013, Heat-Inactivated

Catalog No. NR-51666

Product Description: NR-51666 is a preparation of heat inactivated *Yersinia pestis (Y. pestis)*, strain PY-013 that has been inactivated by heating to 95°C for 20 minutes.

Lot¹: 70025950 Manufacturing Date: 04MAR2019

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology		
Pre-inactivation Post-inactivation	Gram-negative rods No viable bacteria detected	Gram-negative rods No viable bacteria detected
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs) Pre-inactivation Post-inactivation	≥ 99% sequence identity to <i>Y. pestis</i> type strain NCTC 5923 (GenBank: NR_025160.1) ≥ 99% sequence identity to <i>Y. pestis</i> , strain PY-013 (NR-31163)	100% sequence identity to <i>Y. pestis</i> type strain NCTC 5923 (GenBank: NR_025160.1) 100% sequence identity to <i>Y. pestis</i> , strain PY-013 (NR-31163)
Functional Activity by PCR Pre-inactivation Post-inactivation	700 to 1500 base pair amplicon 700 to 1500 base pair amplicon	~ 1500 base pair amplicon ~ 1500 base pair amplicon
SDS-PAGE Analysis (Coomassie Blue Staining)	Protein profile similar to active Y. pestis	Protein profile similar to active Y. pestis (Figure 1) ²
Functional Activity by Western Blot Analysis Antibody to YopM protein, clone CE 15-1 Antibody to YopM protein, clone 2A3.3A8.1A2 Antibody to LcrV protein Antiserum to F-1 antigen	Reactive Reactive Reactive Reactive	Reactive (Figure 2) ^{2,3} Reactive (Figure 3) ^{2,4} Reactive (Figure 4) ^{2,5} Reactive (Figure 5) ^{2,6}
Purity (post-freeze) ⁷	No growth	No growth
Bacterial Inactivation (14 days) 10% of heat-inactivated preparation plated on agar ⁸	No viable bacteria detected	No viable bacteria detected

¹NR-51666 was inactivated by heating *Y. pestis*, strain PY-013 to 95°C for 20 minutes. The material was cooled to room temperature, vialed and frozen to produce this lot. The material used for inactivation was produced from the second passage of BEI Resources NR-31163 lot 61650866, grown under propagation conditions.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

²Results comparable to Y. pestis, strain KIM Derivative 19 (D19) (BEI Resources NR-4681 lot 70012951)

³Using a 1:500 dilution of monoclonal anti-Y. *pestis* outer protein M (YopM), clone 2A3.3A8.1A2 (BEI Resources NR-799) as a primary antibody and a 1:2000 dilution of goat anti-mouse IgG with HRP conjugate as a secondary antibody

⁴Using a 1:500 dilution of monoclonal anti-Y. *pestis* outer protein M (YopM), clone CE 15-1 (BEI Resources NR-800) as a primary antibody and a 1:2000 dilution of goat anti-mouse IgG with HRP conjugate as a secondary antibody

⁵Using a 1:500 dilution of monoclonal anti-Y. *pestis* LcrV protein (BEI Resources NR-3831) as a primary antibody and a 1:2000 dilution of goat anti-mouse IgG with HRP conjugate as a secondary antibody

⁶Using a 1:500 dilution of polyclonal anti-Y. pestis F1-antigen (BEI Resources NR-31024) as a primary antibody and a 1:2000 dilution of rabbit antigoat IgG with HRP conjugate as a secondary antibody

⁷Purity of this lot was assessed for 14 days at 37°C under propagation conditions.


⁸NR-51666 was incubated 14 days under propagation conditions. The inactivation procedure that was used has been shown to consistently inactivate 100% of *Yersinia* spp.



SUPPORTING INFECTIOUS DISEASE RESEARCH

Certificate of Analysis for NR-51666


Figure 1: SDS-PAGE Analysis of Heat-Inactivated Y. pestis, strain PY-013 (NR-51666)



Lane 1: SeeBlue™ Plus2 Pre-Stained Protein Standard

(Invitrogen™; 6 μL) Lane 2: NR-4681 (26 μL) Lane 3: NR-51666 (26 μL)


Figure 2: Western Blot with Y. pestis YopM mAb (NR-799)



Lane 1: SeeBlue™ Plus2 Pre-Stained Protein Standard (6 µL)

Lane 2: NR-4681 (26 μL) Lane 3: NR-51666 (26 μL)


Figure 4: Western Blot with Y. pestis LcrV mAb (NR-3831)



Lane 1: SeeBlue™ Plus2 Pre-Stained Protein Standard (6 µL)

Lane 2: NR-4681 (26 μL) Lane 3: NR-51666 (26 μL)


Figure 3: Western Blot with Y. pestis YopM mAb (NR-800)



Lane 1: SeeBlue™ Plus2 Pre-Stained Protein Standard (6 µL)

Lane 2: NR-4681 (26 μL) Lane 3: NR-51666 (26 μL)

Figure 5: Western Blot with Y. pestis F-1 mAb (NR-31024)



Lane 1: SeeBlue™ Plus2 Pre-Stained Protein Standard (6 µL)

Lane 2: NR-4681 (26 μL) Lane 3: NR-51666 (26 μL)

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898



Certificate of Analysis for NR-51666

/Heather Couch/ <u>Heather Couch</u>

16 JUL 2019

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

NR-51666_70025950_16JUL2019