

***Mycobacterium tuberculosis* Gateway® Clone Set, Recombinant in *Escherichia coli*, Plate 10**

Catalog No. NR-19646

This reagent is the tangible property of the U.S. Government.

For research use only. Not for human use.

Contributor:

Pathogen Functional Genomics Resource Center at the J. Craig Venter Institute

Manufacturer:

BEI Resources

Product Description:

Clone plates are replicated using a BioMek® FX robot. Production in the 96-well format has increased risk of cross-contamination between adjacent wells. Individual clones should be purified (e.g. single colony isolation and purification using good microbiological practices) and sequence-verified prior to use. BEI Resources only confirms the clone plate orientation and viability of randomly picked clones. BEI Resources does not confirm or validate individual clone identities provided by the contributor.

The *Mycobacterium tuberculosis* (*M. tuberculosis*) Gateway® clone set consists of 42 plates which contain 3724 sequence validated clones (3294 *M. tuberculosis*, strain H37Rv clones supplemented with 430 unique open reading frames (ORF) from *M. tuberculosis*, strain CDC1551) cloned in *Escherichia coli* (*E. coli*) DH10B-T1 cells. Each ORF was recombined in vector pDONR™221 with an ATG start codon and no stop codon. The sequence was validated by full length sequencing of each entry clone with greater than 1X coverage and a mutation rate of less than 0.2%. Detailed information about each clone is shown in Table 1.

Information related to the use of Gateway® Clones can be obtained from [Invitrogen™](#). Recombination was facilitated through an *attB* substrate (*attB*-PCR product or a linearized *attB* expression clone) with an *attP* substrate (pDONR™221) to create an *attL*-containing entry clone. The entry clone contains recombinational cloning sites, *attL1* and *attL2* to facilitate gene transfer into a destination vector, M13 forward and reverse priming sites for sequencing and a kanamycin resistance gene for selection. Please refer to the [Invitrogen™ Gateway® Technology Manual](#) for additional details.

Plate orientation and viability were confirmed for NR-19646.

Material Provided:

Each inoculated well of the 96-well plate contains approximately 60 µL of *E. coli* culture (strain DH10B-T1) in Luria Bertani (LB) Broth containing 50 µg/mL kanamycin supplemented with 15% glycerol.

Packaging/Storage:

NR-19646 was packaged aseptically in a 96-well plate. The product is provided frozen and should be stored at -80°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

Growth Conditions:

Media:

LB Broth or Agar containing 50 µg/mL kanamycin.

Incubation:

Temperature: *E. coli*, strain DH10B-T1 clones should be grown at 37°C.

Atmosphere: Aerobic

Propagation:

1. Scrape top of frozen well with a pipette tip and streak onto agar plate.
2. Incubate the plates at 37°C for 18 to 24 hours.

Citation:

Acknowledgment for publications should read “The following reagent was obtained through BEI Resources, NIAID, NIH: *Mycobacterium tuberculosis* Gateway® Clone Set, Recombinant in *Escherichia coli*, Plate 10, NR-19646.”

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. [Biosafety in Microbiological and Biomedical Laboratories](#). 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see www.cdc.gov/biosafety/publications/bmb15/index.htm.

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI

Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

References:

1. Cole, S. T., et al. "Deciphering the Biology of *Mycobacterium tuberculosis* from the Complete Genome Sequence." *Nature* 393 (1998): 537-544. PubMed: 9634230.
2. Camus, J. C., et al. "Re-Annotation of the Genome Sequence of *Mycobacterium tuberculosis* H37Rv." *Microbiology* 148 (2002): 2967-2973. PubMed 12368430.

ATCC® is a trademark of the American Type Culture Collection



Table 1: *Mycobacterium tuberculosis*, Gateway® Clones, Plate 10 (ZMTDJ)¹

Clone	Well Position	ORF Length	Locus ID	Description	Accession Number	Average Depth of Coverage
76977	A01	664	MT1082	hypothetical protein	NP_335520.1	1.605421687
77029	A02	673	MT0903	hypothetical protein	NP_335331.1	1.930163447
77104	A03	673	MT3084	hypothetical protein	NP_337599.1	-
77020	A04	679	MT3503	hypothetical protein	NP_338026.1	2
76973	A05	685	MT0726	hypothetical protein	NP_335142.1	2
77056	A06	715	MT0177	hypothetical protein	NP_334583.1	2
77034	A07	748	MT2934	hypothetical protein	NP_337445.1	3.965240642
77150	A08	748	MT3770	hypothetical protein	NP_338323.1	3.881016043
77001	A09	847	MT3393	hypothetical protein	NP_337922.1	1.551357733

¹All information in this table was provided by J. Craig Venter Institute at the time of deposition.