**b**|**e**|**i** resources

SUPPORTING INFECTIOUS DISEASE RESEARCH

# *Streptococcus pneumoniae* Gateway<sup>®</sup> Clone Set, Recombinant in *Escherichia coli*, Plate 16

# Catalog No. NR-19583

This reagent is the tangible property of the U.S. Government.

# For research use only. Not for human use.

#### **Contributor:**

Pathogen Functional Genomics Resource Center at the J. Craig Venter Institute

#### Manufacturer:

**BEI Resources** 

#### **Product Description:**

Clone plates are replicated using a BioMek<sup>®</sup> FX robot. Production in the 96-well format has increased risk of crosscontamination between adjacent wells. Individual clones should be purified (e.g. single colony isolation and purification using good microbiological practices) and sequence-verified prior to use. BEI Resources only confirms the clone plate orientation and viability of randomly picked clones. BEI Resources does not confirm or validate individual clone identities provided by the contributor.

The Streptococcus pneumoniae (S. pneumoniae) Gateway<sup>®</sup> clone set consists of approximately 2029 sequence validated clones from S. pneumonia, strain TIGR4 cloned in *Escherichia coli* (*E. coli*) DH10B-T1 cells.<sup>1</sup> Each open reading frame was constructed in vector <u>pDONR<sup>TM</sup>221</u> (Invitrogen<sup>TM</sup>) with a native start codon and no stop codon. The sequence was validated by full length sequencing of each clone with greater than 1X coverage and a mutation rate of less than 0.2%. Detailed information about each clone is shown in Table 1.

Information related to the use of Gateway<sup>®</sup> Clones can be obtained from Invitrogen<sup>TM</sup>. Recombination was facilitated through an *att*B substrate (*att*B-PCR product or a linearized *att*B expression clone) with an *att*P substrate (pDONR<sup>TM</sup>221) to create an *att*L-containing entry clone. The entry clone contains recombinational cloning sites, *att*L1 and *att*L2 to facilitate gene transfer into a destination vector, M13 forward and reverse priming sites for sequencing and a kanamycin resistance gene for selection. Please refer to the Invitrogen<sup>TM</sup> <u>Gateway<sup>®</sup> Technology Manual</u> for additional details.

Plate orientation and viability were confirmed for NR-19583.

# Material Provided:

Each inoculated well of the 96-well plate contains approximately 60  $\mu$ L of *E. coli* culture (strain DH10B-T1) in Luria Bertani (LB) broth containing 50  $\mu$ g/mL kanamycin

supplemented with 15% glycerol.

# Packaging/Storage:

NR-19583 was packaged aseptically in a 96-well plate. The product is provided frozen and should be stored at -80°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

# **Growth Conditions:**

### Media:

LB broth containing 50 μg/mL kanamycin LB agar containing 50 μg/mL kanamycin <u>Incubation</u>: Temperature: 37°C. Atmosphere: Aerobic <u>Propagation</u>: 1. Scrape top of frozen well with a pipette tip and streak onto agar plate.

2. Incubate the plates at 37°C for 18 to 24 hours.

### Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: *Streptococcus pneumoniae* Gateway<sup>®</sup> Clone Set, Recombinant in *Escherichia coli*, Plate 16, NR-19583."

# **Biosafety Level: 1**

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. <u>Biosafety in</u> <u>Microbiological and Biomedical Laboratories</u>. 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see www.cdc.gov/biosafety/publications/bmbl5/index.htm.

#### **Disclaimers:**

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at <u>www.beiresources.org</u>.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC<sup>®</sup> nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC<sup>®</sup> nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal.  $ATCC^{\otimes}$  and the U.S. Government are not liable for any

E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898 **b**|**e**|**i** resources

SUPPORTING INFECTIOUS DISEASE RESEARCH

damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure Government, ATCC<sup>®</sup>, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

#### **Use Restrictions:**

This material is distributed for internal research, noncommercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale. authenticity and reliability of materials on deposit, the U.S.

### **References:**

 Kwon, K., et al. "A Correlation Analysis of Protein Characteristics Associated with Genome-Wide High Throughput Expression and Solubility of *Streptococcus pneumoniae* Proteins." <u>Protein Expr. Purif.</u> 55 (2007): 368-378. PubMed: 17703947.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection.



#### Table 1: Streptococcus pneumoniae Gateway<sup>®</sup> Clone Set, Recombinant in Escherichia coli, Plate 16 (YSPCQ)<sup>1</sup>

Clone	Well	Locus ID	Description	ORF	Accession	Average Depth of
67521		SD1335	hypothetical protein SP 1335	1/10	ND 345703 1	20.825
67420	A01	SP2019	transposaso	1410	NF_040790.1	20.023
67245	A02	SP2010	hypothetical protain SD 0011	970	- ND 245206 1	2.0152
67202	A03	SF0911	hypothetical protein SF_0911	070	NP_343390.1	0.0931
07393	A04	SP0606	hypothetical protein SP_0866	037	NP_345353.1	0.1004
67509	A05	SP0068	h setaria sin Dial	195	NP_344617.1	3.8923
67525	A06	SP0531		561	NP_345048.1	3.4866
67317	A07	SP1514	FUF1 ATP synthase subunit C	193	NP_345965.1	3
67281	A08	SP0134	hypothetical protein SP_0134	243	NP_344677.1	12.988
67421	A09	SP1904	hypothetical protein SP_1904	532	NP_346335.1	2
67221	A10	SP0584	transcriptional regulator	609	NP_345097.1	5.0279
67438	A11	SP0451	hypothetical protein SP_0451	504	NP_344972.1	3.1171
67397	A12	SP2119	transcriptional regulator	1528	NP_346538.1	3.411
67873	B01	SP1065	hypothetical protein SP_1065	1020	NP_345539.1	10.931
67569	B02	SP0418	acyl carrier protein	477	NP_344941.1	4.3375
67577	B03	SP1948	hypothetical protein SP_1948	766	NP_346376.1	2
67905	B04	SP0055	hypothetical protein SP_0055	180	NP_344604.1	3.95
67790	B05	SP1694	hypothetical protein SP_1694	919	NP_346132.1	4.3079
67766	B06	SP0041	bacteriocin BlpU	114	NP_344590.1	-
67865	B07	SP0840	hypothetical protein SP_0840	810	NP_345331.1	9.1074
67729	B08	SP1594	hypothetical protein SP_1592	406		4
67693	B09	SP0116	hypothetical protein SP_0116	235	NP_344662.1	14.936
67845	B11	SP0539	bacteriocin BlpM	570	NP_345056.1	4.7175
67621	B12	SP0238	hypothetical protein SP_0238	318	NP_344778.1	5.6352
67534	C01	SP0733	hypothetical protein SP_0733	711	NP_345233.1	12.105
67797	C02	SP1113	DNA-binding protein HU	1095	NP_345584.1	10.742
67637	C03	SP1349	hypothetical protein SP_1349	1473	NP_345807.1	10.688
67749	C04	SP1404	hypothetical protein SP_1404	1800	NP_345862.1	13.55
67841	C05	SP1048	hypothetical protein SP 1048	1005	NP 345522.1	10.816
67753	C06	SP0471	hypothetical protein SP_0471	519	NP_344990.1	2.9114
67586	C07	SP0109	bacteriocin	231	NP_344656.1	15
67785	C08	SP2029	preprotein translocase subunit YaiC	1069		1.8447
67757	C09	SP0208	30S ribosomal protein S10	300	 NP_344748.1	4.62
67821	C10	SP0492	hypothetical protein SP_0492	537	NP_345010.1	4.0205

**b**|**e**|**i** resources

# **Product Information Sheet for NR-19583**

SUPPORTING INFECTIOUS DISEASE RESEARCH

Clone	Well	Locus ID	Description	ORF	Accession	Average Depth of
67901		SD0910	IS620 Spp1_trapapagag Orf2	<b>Lengin</b>		0 1126
67561	C12	SP0019		1020	NF_343311.1	9.1130
67640	D01	SF 1000	ISTS01, transposase OTA	1020	- ND 245759 1	0.0200
67049	D01	SP 1294	by nother tipel protein CD 0882	040	NP_343736.1	12.440
07037	D02	SP0663	hypothetical protein SP_0003	049	NP_343370.1	7.0090
67597	D03	SP0207	hypothetical protein SP_0207	298	NP_344747.1	5.0840
07089	D04	SP1487	hypothetical protein SP_1487	163	NP_345939.1	4
67665	D05	SP1708	hypothetical protein SP_1708	970	NP_346146.1	4.3711
67745	D06	SP0990	nypotnetical protein SP_0990	948	NP_345469.1	9.1097
67593	D07	SP1744	Iojap-related protein	1297	NP_346180.1	3.3192
67629	D08	SP0412	hypothetical protein SP_0412	4/4	<u>NP_344935.1</u>	4.3418
67901	D09	SP1218	nypothetical protein SP_1218	1227	NP_345684.1	7.6365
67725	D10	SP2042	ribonuclease P	1192	NP_346467.1	2.9069
67793	D11	SP0962	lactoylglutathione lyase	921	NP_345443.1	6.9989
67869	D12	SP1334	hypothetical protein SP_1334	1407	NP_345792.1	20.1
67733	E01	SP1567	endoribonuclease L-PSP	283	NP_346014.1	4
68037	E02	SP0941	IS1381 transposase protein B	897	NP_345424.1	9.0011
68069	E03	SP0968	diacylglycerol kinase	927	NP_345449.1	8.5566
68074	E04	SP2074	transposase	1312	-	1.3323
68053	E05	SP0189	hypothetical protein SP_0188	288	NP_344730.1	4.2326
67945	E06	SP1310	IS1381 transposase protein A	1344	NP_345771.1	11.839
68249	E07	SP1851	hypothetical protein SP_1851	349	NP_346283.1	2
68005	E08	SP0015	IS630-Spn1, transposase Orf1	102	NP_344568.1	-
68161	E09	SP2015	IS630-Spn1, transposase Orf1	988	NP_346442.1	1.9706
67917	E10	SP0424	(3R)-hydroxymyristoyl-ACP dehydratase	480	NP_344947.1	4.2583
68241	E11	SP0433	transcription antitermination protein NusB	489	NP_344955.1	5.7546
68050	E12	SP0817	MutT/nudix family protein	792	NP_345310.1	9.5859
68257	F01	SP1140	hypothetical protein SP_1140	1137	NP_345610.1	9.6403
68145	F02	SP0416	MarR family transcriptional regulator	474	NP_344939.1	2.981
68265	F03	SP1636	Rrf2 family protein	664	NP_346076.1	5.6114
68030	F04	SP0849	-	822	-	8.9197
68117	F05	SP0229	50S ribosomal protein L15	312	NP_344769.1	5.25
68129	F06	SP0135	glycosyl transferase	246	NP_344678.1	11.963
68217	F07	SP1297	flavodoxin	1338	NP_345761.1	10.918
68237	F08	SP2165	fucose operon FucU protein	1861	NP_346579.1	3.5943
67981	F09	SP1127	hypothetical protein SP_1127	1116	NP_345597.1	10.571
68102	F10	SP1520	acetyltransferase	196	NP_345970.1	3
67929	F11	SP0992	hypothetical protein SP_0992	951	NP_345471.1	8.4101
68285	F12	SP0024	hypothetical protein SP_0024	105	NP_344575.1	2
68137	G01	SP0811	transposase family protein	786	NP_345307.1	9.374
68090	G02	SP0507	type I restriction-modification system, S subunit	546	NP_345025.1	7.0678
68645	G03	SP0719	hypothetical protein SP_0719	702	NP_345220.1	13.729
68473	G04	SP1773	-	142	-	1.9014
68457	G05	SP1346	hypothetical protein SP 1346	1455	NP_345804.1	15.96
68489	G06	SP0205	anaerobic ribonucleoside-triphosphate reductase activating protein	297	 NP_344745.1	4.2323
68593	G07	SP0179	Holliday junction DNA helicase RuvA	282	NP 344720.1	3.1596
68357	G08	SP0370	Holliday junction-specific endonuclease	444	NP_344897.1	3.3153

bei resources

# **Product Information Sheet for NR-19583**

SUPPORTING INFECTIOUS DISEASE RESEARCH

Clone	Well Position	Locus ID	Description	ORF Length	Accession Number	Average Depth of Coverage
68369	G09	SP1170	hypothetical protein SP_1170	1176	NP_345639.1	8.881
68477	G10	SP1041	hypothetical protein SP_1041	996	NP_345516.1	7.6596
68401	G11	SP1625	cadmium resistance transporter	621	NP_346065.1	3.0628
68617	G12	SP0589	serine acetyltransferase	609	NP_345102.1	3.197
68553	H01	SP1134	hypothetical protein SP_1134	1128	NP_345604.1	10.676
68493	H02	SP1007	hypothetical protein SP_1007	960	NP_345483.1	8.974
68521	H03	SP1601	hypothetical protein SP_1601	433	NP_346045.1	3.8891
68501	H04	SP0104	HAD superfamily hydrolase	228	NP_344651.1	17.943
68373	H05	SP0725	thiamine-phosphate pyrophosphorylase	705	NP_345226.1	13.691
68301	H06	SP1957	ABC transporter ATP-binding protein	793	NP_346384.1	2
68417	H07	SP0695	HesA/MoeB/ThiF family protein	690	NP_345200.1	11.038
68333	H08	SP0687	ABC transporter ATP-binding protein	681	NP_345192.1	13.088
68565	H09	SP0851	glycerol-3-phosphate acyltransferase PIsY	825	NP_345340.1	8.6618
68537	H10	SP0489	PAP2 family protein	534	NP_345007.1	2.191
68413	H11	SP0721	hypothetical protein SP_0721	702	NP_345222.1	12
68601	H12	SP1992	cell wall surface anchor family protein	883	NP_346419.1	2

<sup>1</sup>All information in this table was provided by J. Craig Venter Institute at the time of deposition.