

Certificate of Analysis for HM-687

Bacteriophage Phi05 1315B/Pa1400, Infectious for Pseudomonas aeruginosa

Catalog No. HM-687

Product Description: Bacteriophage phi05_1315B/Pa1400 was derived from a blood sample taken from a patient with septicemia. Bacteriophage phi05_1315B/Pa1400 (HMP ID 9798) is a reference genome for The Human Microbiome Project (HMP).

Lot^{1,2}: 60277714 Manufacturing Date: 26AUG2011

| TEST | SPECIFICATIONS | RESULTS |
|---|-----------------------------|--------------------------------|
| Plaque Morphology ³ | Report results | Small and clear |
| Titer ³ Plaque-forming units (pfu) with Pseudomonas aeruginosa | >10 ⁶ pfu per mL | 2 x 10 ⁹ pfu per mL |
| Sterility | 0.22 µm filtered | 0.22 µm filtered |

¹Quality control of HMP material is only performed to demonstrate that the product distributed by BEI Resources is identical to the deposited material. It should not be considered a complete characterization of the deposited material.

Date: 17 NOV 2014

Signature: (

Title:

Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²The deposited material was added to a Luria-Burtani (LB) soft agar overlay (0.5%) supplemented with 10 mM MgSO₄ containing the phage host, *Pseudomonas aeruginosa*, strain Pa1400. The soft agar overlay was added to an LB agar kolle supplemented with 10 mM MgSO₄ and incubated 24 hours at 30°C in an aerobic atmosphere. After one additional passage under the above conditions, the bacteriophage were harvested to produce this lot.

³Lysis was observed on the soft agar overlay after 24 hours at 30°C in an aerobic atmosphere.