

SUPPORTING INFECTIOUS DISEASE RESEARCH

Mycobacterium tuberculosis, Strain H37Rv, Gamma-Irradiated Whole Cells

Catalog No. NR-14819

This reagent is the tangible property of the U.S. Government.

Product Description: *Mycobacterium tuberculosis*, strain H37Rv, was grown to late-log phase in glycerol-alanine-salts medium and inactivated by exposure to 2.4 megaRads of ionizing gamma irradiation using a ¹³⁷Cs source.

Lot: 61757227

Manufacturing Date: 26AUG2014

Production and QC testing were performed by Colorado State University. The Colorado State University documentation for bulk lot 14.Rv.2.7.22.7.WCg.a is attached.

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by the contractor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.



WORK SHEET FOR γ-IRRADIATED CELLS

General Information

Lot Number: <u>14.Rv.2.7.22.7.WCg.a</u> Species: M. tuberculosis Strain: H37Rv

Growth Conditions: GAS Medium:___ # of Fernbachs: <u>8</u> # of Roller bottles: 40 Inoculation Date: 7.22.14 Harvest Date: 8.5.14 SOP # used for Harvest: PP003.3 1/11/13 version

Wet Weight: 96.2g

Date irradiated: <u>8.26.14</u> By: JB

Alamar Blue QC Information

Date Started: 9.17.14 Date Checked: 9.19.14

Positive Control: Positive X Negative

Sample: Positive Negative X

Aliquot Information

96.2g in 230ml falcon tube 1 x

10-31-14 (date) (Research Associate)

(Laboratory Supervisor)

<u>10/3</u>//04 (date)

| | | Initial 0 hr 570/600 | 570/600 | | |
|---------|---------------------|----------------------|---------------------|--------|---------|
| | 9-17-14. | 0.77 | 0.74 | 0.70 | 0.72 |
| N.2.7.5 | 4.Rv.2.7.22.7.WCg.a | 0.78 | 0.76 | 0.72 | 0.72 |
| | | 0.78 | 0.76 | 0.75 | 0.73 |
| | | 0.77 | 0.79 | 0.76 | 0.74 |
| | | | | | |
| | | Initia | Initial Minus Blank | × | |
| | | 0D.1 | 0D .01 | 0D.001 | OD.0001 |
| | | 0.04 | 0.01 | -0.03 | -0.01 |
| | Sample | 0.06 | • 0.04 | 0.00 | 00.00 |
| | | 0.06 | 0.04 | 0.02 | 00.00 |
| | Pos. Cont. | 0.04 | 0.06 | 0.03 | 0.01 |
| | | | | | |
| | | ini | Initial Average | | |
| | | 0D.1 | 0D .01 | 0D.001 | OD.0001 |
| | Sample | 0.05 | 0.03 | 0.00 | 00.0 |
| e | Pos. Ctrl. | 0.04 | 0.06 | 0.03 | 0.01 |

| | Idual | e .01 | Sample .01 Pos .1 | Pos .01 |
|---|-------|-------|-------------------|---------|
| 0 | 0.05 | 0.03 | 0.04 | 0.06 |
| 0 | 0.08 | 0.01 | 2.79 | 0.67 |
| - | 0.13 | 0.02 | 3.30 | 1.21 |

| | | 0.74 | 0.74 | 0.76 | | OD.0001 | |
|--------------|----------|------|------|------|------------------|---------|--|
| | 0.73 | 0.74 | 0.74 | 0.80 | × | 0D .001 | |
| /600 | 0.74 | 0.74 | 0.75 | 1.40 | 24hr Minus Blank | 0D .01 | |
| 24hr 570/600 | 0.80 | 0.81 | 0.82 | 3.53 | 24hr | 0D.1 | |
| | 9-18-14. | | | | | | |
| L | | | | | | | |

| | 0.07 | 0.01 | 00.00 | 0.00 |
|------------|------|-------|---------|---------|
| Sample | 0.08 | 0.01 | 0.00 | 0.00 |
| | 0.09 | 0.01 | 0.00 | 0.00 |
| Pos. Cont. | 2.79 | 0.67 | 0.07 | 0.02 |
| | 0D.1 | 0D.01 | 0D .001 | OD.0001 |
| Sample | 0.08 | 0.01 | 00.0 | 0.00 |
| Pos. Ctrl. | 2.79 | 0.67 | 0.07 | 0.02 |

| | 000/0/C 111 24 | PUU | | |
|----------|----------------|------|------|------|
| 9-19-14. | 0.85 | 0.76 | 0.74 | 0.74 |
| | 0.87 | 0.76 | 0.75 | 0.75 |
| | 06.0 | 0.76 | 0.75 | 0.75 |
| | 4.05 | 1.95 | 0.84 | 0.78 |

| | 0D.1 | 0D .01 | OD .001 OD .0001 | OD.0001 |
|-----------|------|--------|------------------|---------|
| | 0.11 | 0.01 | 00.00 | 0.00 |
| Sample | 0.13 | 0.02 | 0.01 | 0.00 |
| | 0.15 | 0.02 | 0.01 | 0.00 |
| os. Cont. | 3.30 | 1.21 | 0.10 | 0.03 |

| | 48 | 48hr Average | | |
|------------|------|--------------|---------|---------|
| | 0D.1 | 0D .01 | 0D .001 | OD.0001 |
| Sample | 0.13 | 0.02 | 0.00 | 0.00 |
| Pos. Ctrl. | 3.30 | 1.21 | 0.10 | 0.03 |

