

## *Naegleria fowleri*, Strain CDC:V413

### Catalog No. NR-46493

This reagent is the tangible property of the U.S. Government.

### For research use only. Not for human use.

#### Contributor:

Govinda S. Visvesvara, Ph.D., and Michael Arrowood, Ph.D., Centers for Disease Control and Prevention, National Center for Emerging and Zoonotic Infectious Diseases, Division of Foodborne, Waterborne and Environmental Diseases, Waterborne Disease Prevention Branch, Atlanta, Georgia, USA

#### Manufacturer:

BEI Resources

#### Product Description:

Protozoa Classification: *Vahlkampfiidae*, *Naegleria*

Species: *Naegleria fowleri*

Strain: CDC:V413

Original Source: *Naegleria fowleri* (*N. fowleri*), strain

CDC:V413 is a clinical isolate collected in 1998 from the cerebrospinal fluid of a 17-year-old man in Texas, USA.<sup>1,2</sup>

Comment: *N. fowleri*, strain CDC:V413 was deposited to BEI Resources as genotype I.<sup>1,2</sup>

*N. fowleri* is a free-living pathogenic amoeboflagellate that feeds mainly on bacteria, yeast and algae. It is the causative agent of primary amoebic meningoencephalitis, a water-borne disease of the central nervous system in humans and is associated with recreational activities in contaminated waters.<sup>2-8</sup> Infection occurs when the amoeba enters the human nasal cavity where it attaches to the nasal mucosa and travels along olfactory nerves, eventually entering the brain, where it causes extensive tissue damage, inflammation, and hemorrhagic necrosis.<sup>5</sup> *N. fowleri* is a moderate thermophile and can tolerate temperatures up to 45°C. It has been isolated worldwide from soil and fresh waters naturally heated by the sun, including lakes, ponds, well water, geothermal springs and in areas thermally polluted by industries.<sup>2-8</sup>

The *N. fowleri* life cycle consists of three morphological stages: a dividing, feeding, infective amoeboid trophozoite; a transitory pear-shaped flagellate formed from the amoeba during conditions of nutrient deprivation in water; and, a resistant cyst that forms during adverse environmental conditions such as food deprivation, crowding, ionic changes and the presence of toxin-producing bacteria.<sup>5,6,8</sup> Eight genotypes of *N. fowleri* have been identified using the internal transcribed spacer regions and 5.8S ribosomal RNA gene sequences.<sup>2-5</sup>

#### Material Provided:

Each vial of NR-46493 contains approximately 0.5 mL of culture in cryopreservative [7.5% dimethylsulfoxide (DMSO)]. Please refer to the Certificate of Analysis for the specific

culture media used for each lot and refer to Appendix I for cryopreservation instructions.

#### Packaging/Storage:

NR-46493 was packaged aseptically in screw-capped plastic cryovials and is provided frozen on dry ice. The product should be stored at cryogenic temperature (-130°C or colder), preferably in the vapor phase of a liquid nitrogen freezer. If liquid nitrogen storage facilities are not available, frozen cryovials may be stored at -70°C or colder for approximately one week.

Note: Do not under any circumstances store vials at temperatures warmer than -70°C. Storage under these conditions will result in the death of the culture.

To insure the highest level of viability, the culture should be initiated immediately upon receipt. Any warming of the product during shipping and transfer must be avoided, as this will adversely affect the viability of the product. For transfer between freezers and for shipping, the product may be placed on dry ice for brief periods, although use of a portable liquid nitrogen carrier is preferred. Please read the following recommendations prior to using this material.

#### Growth Conditions:

Modified Peptone - Yeast Extract - Nucleic Acid - Folic Acid - Hemin (PYNFH) medium (ATCC® Medium 1034) supplemented with 10% heat-inactivated fetal bovine serum (Appendix II)

#### Incubation:

Temperature: 35°C

Atmosphere: Aerobic

#### Propagation:

1. Place the frozen vial in a 35°C to 37°C water bath and thaw for approximately 2 to 3 minutes. Do not agitate the vial. Do not leave the vial in the water bath after it is thawed.
2. Immediately after thawing, aseptically transfer the contents of the vial to a T-25 tissue culture flask containing 5 to 10 mL modified PYNFH medium.
3. Screw the cap on tightly and incubate the tube or flask at 35°C.

#### Maintenance:

1. When the culture is at or near peak density, vigorously agitate or scrape the surface of the flask to detach adherent cells.
2. Transfer approximately 0.25 mL to a new flask containing 5 to 10 mL of freshly made modified PYNFH medium.
3. Screw the caps on tightly and incubate at 35°C.
4. The amoeba will form an almost continuous sheet of cells on the bottom surface of the flask. Repeat steps 1 through 3 every 7 to 10 days.

Please refer to Appendix I for cryopreservation instructions.

#### Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: *Naegleria fowleri*, Strain CDC:V413, NR-46493."

**Biosafety Level: 2**

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see [www.cdc.gov/biosafety/publications/bmbl5/index.htm](http://www.cdc.gov/biosafety/publications/bmbl5/index.htm).

**Disclaimers:**

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at [www.beiresources.org](http://www.beiresources.org).

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

**Use Restrictions:**

**This material is distributed for internal research, non-commercial purposes only.** This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

**References:**

1. Visvesvara, G. S., Personal Communication.
2. Zhou, L., et al. "Genetic Variations in the Internal Transcribed Spacer and Mitochondrial Small Subunit rRNA Gene of *Naegleria* Spp." J. Eukaryot. Microbiol. 50 (2003): 522-526. PubMed: 14736150.
3. De Jonckheere, J. F. "Origin and Evolution of the Worldwide Distributed Pathogenic Amoeboflagellate *Naegleria fowleri*." Infect. Genet. Evol. 11 (2011): 1520-1528. PubMed: 21843657.
4. Marciano-Cabral, F. and G. A. Cabral. "The Immune Response to *Naegleria fowleri* Amebae and Pathogenesis

- of Infection." FEMS Immunol. Med. Microbiol. 51 (2007): 243-259. PubMed: 17894804.
5. Pélandakis M., et al. "Analysis of the 5.8S rRNA Gene and the Internal Transcribed Spacers in *Naegleria* Spp. and in *N. fowleri*." J. Eukaryot. Microbiol. 47 (2000): 116-121. PubMed: 10750838.
6. Marciano-Cabral, F. "Biology of *Naegleria* Spp." Microbiol. Rev. 52 (1988): 114-133. PubMed: 3280964.
7. De Jonckheere, J. F. "A Century of Research on the Amoeboflagellate Genus *Naegleria*." Acta Protozool. 41 (2002): 309-342.
8. Visvesvara, G. S., H. Moura and F. L. Schuster. "Pathogenic and Opportunistic Free-Living Amoebae: *Acanthamoeba* Spp., *Balamuthia mandrillaris*, *Naegleria fowleri*, and *Sappinia dipoloidea*." FEMS Immunol. Med. Microbiol. 50 (2007): 1-26. PubMed: 17428307.
9. Martinez-Castillo, M., et al. "*Naegleria fowleri* after 50 Years: Is it a Neglected Pathogen?" J. Med. Microbiol. 65 (2016): 885-896. PubMed: 27381464.

ATCC® is a trademark of the American Type Culture Collection.



**APPENDIX I: CRYOPRESERVATION**

1. Harvest *Naegleria* from multiple flasks by scraping the surface of the flask with a cell scraper to detach adhering trophozoites.
2. Transfer the cell suspensions to 15 mL or 50 mL plastic centrifuge tubes.
3. Adjust the cell concentration to  $1 \times 10^6$  to  $2 \times 10^7$  cells/mL with fresh modified PYNFH medium.  
Note: If the concentration of cells is too low, centrifuge at  $1300 \times g$  for 10 minutes and resuspend in a smaller volume of fresh medium to yield the desired cell concentration.
4. Mix equal volumes of cell suspension and fresh medium containing 15% DMSO to yield a final concentration of  $1 \times 10^6$  to  $2 \times 10^7$  cells/mL in 7.5% DMSO. The freezing process should start 15 to 30 minutes following the addition of cryoprotective solution to the cell suspension.
5. Dispense 0.5 mL aliquots into 1 to 2 mL sterile plastic screw-capped vials for cryopreservation.
6. Place the vials in a controlled rate freezing unit. From room temperature cool the vials at  $-1^\circ\text{C}/\text{min}$  to  $-40^\circ\text{C}$ . If the freezing unit can compensate for the heat of fusion, maintain rate at  $-1^\circ\text{C}/\text{min}$  through this phase. At  $-40^\circ\text{C}$ , plunge vials into liquid nitrogen. Alternatively, place the vials in a Nalgene  $1^\circ\text{C}$  freezing container. Place the container at  $-80^\circ\text{C}$  for 1.5 to 2 hours and then plunge vials into liquid nitrogen.
7. Store in either the vapor or liquid phase of a nitrogen refrigerator ( $-130^\circ\text{C}$  or colder).

**APPENDIX II: Modified Peptone - Yeast Extract - Nucleic Acid - Folic Acid - Hemin (PYNFH) Medium (ATCC® MEDIUM 1034)**

1. Prepare the Basal medium (see recipe below), autoclave for 20 minutes at  $121^\circ\text{C}$ , and allow to cool.
2. Prepare the Buffer solution (see recipe below) and filter sterilize.
3. Aseptically prepare the PYNFH medium (see recipe below), mix thoroughly and store at  $4^\circ\text{C}$ :

<u>Basal Medium</u>		<u>Buffer Solution</u>	
Peptone	10.0 g	KH <sub>2</sub> PO <sub>4</sub>	18.1 g
Yeast Extract	10.0 g	Na <sub>2</sub> HPO <sub>4</sub>	25.0 g
Yeast Nucleic Acid	1.0 g	Distilled Water	1.0 L
Folic Acid	15.0 mg		
Hemin	1.0 mg		
Distilled water	880.0 mL		

4. Aseptically prepare the PYNFH medium (see recipe below), mix thoroughly and store at  $4^\circ\text{C}$ :

<u>PYNFH Medium</u>	
Basal medium	880.0 mL
Buffer Solution	20.0 mL

5. Aseptically supplement PYNFH medium with 100.0 mL heat-inactivated fetal bovine serum prior to use.