

## NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

## DATA SHEET

HIV-1 QH0692 Env Expression Vector (QH0692, clone 42 (SVPB6)) Reagent:

Catalog Number: 11018

Lot Number: 050734

**Release Category:** В

Provided: 1 vial containing 20 μg of plasmid DNA in TE buffer (0.5 mg/ml)

**Cloning Vector:** pcDNA3.1D/V5-His TOPO@. The size of the cloning vector including the insert is 8566 bp.

**Cloning Site:** The HIV-1 env/rev cassette was directly cloned into the cloning site of pcDNA3.1D/V5-His

TOPO© expression vector, in the correct orientation with the CMV promoter. The size of

the insert is 3052 bp.

AY835439 GenBank:

**Description:** A PCR fragment containing full-length env and rev genes was derived from the genomic

DNA of infected PBMC. Original virus was obtained by PBMC co-culture. The env/rev cassette was cloned into pcDNA3.1D/V5-His TOPO© expression vector by a directional cloning approach. A single transformed ampicillin resistant E. coli colony was selected and expanded. Recombinant plasmid carries resistance genes for ampicillin and

neomycin. Sequence information is available.

Click here to see the sequence file

Special Characteristics:

Plasmid Expansion: It is recommended that this plasmid be expanded using DH5a™ Competent Cells in LB medium at 34°C. The clone represents env/rev sequences from a subject with acute subtype B infection (female-male transmission in Trinidad/Tobago). The clone expresses a functional env/rev cassette and can be used to generate

pseudotyped infectious virions that use CCR5 as the viral co-receptor. OH0692.42 Env containing pseudovirions are included in a standard virus neutralization panel for subtype

B strains (SVPB6).

This clone is also available as a member of a panel set, see catalog #11227. Click here

to see the table of individual reagents for this reference panel.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 10/18/2017 Page 1 of 2 Recommended Storage:

-80°C.

**Contributor:** Dr. David Montefiori, Dr. Feng Gao and Dr. Ming Li.

**References:** Li, M., Gao F., Mascola J.R., Stamatatos L., Polonis V.R., Koutsoukos M., Voss G.,

Goepfert P., Gilbert P., Greene K.M., Bilska M., Kothe D.L., Salazar-Gonzalez J.F., Wei X., Decker J.M., Hahn B.H., and Montefiori D.C. Human immunodeficiency virus type 1 env clones from acute and early subtype B infections for standardized assessments of vaccine-elicited neutralizing antibodies. *J. Virology* **79**(16): 10108-10125, 2005.

NOTE: Acknowledgment for publications should read "The following reagent was obtained

through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH from Dr. David Montefiori and Dr. Feng Gao: HIV-1 QH0692 Env Expression Vector (QH0692, clone 42

(SVPB6))." Also include the reference cited above in any publications.

**Last Updated:** October 18, 2017

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 10/18/2017 Page 2 of 2