



## NIH AIDS Reagent Program

20301 Century Boulevard  
Building 6, Suite 200  
Germantown, MD 20874  
USA

Phone: 240 686 4740  
Fax: 301 515 4015  
aidsreagent.org

### DATA SHEET

|                                 |   |
|---------------------------------|---|
| <b>Reagent:</b>                 | Human CCR1 Expression Vector (pcCCR1)   |
| <b>Catalog Number:</b>          | 3321  |
| <b>Lot Number:</b>              | 97058   |
| <b>Release Category:</b>        | C   |
| <b>Provided:</b>                | 1 µg of dried purified DNA stabilized in DNastable <i>PLUS</i>  |
| <b>Cloning Site:</b>            | BamHI and XhoI (5'-3'). The size of the insert is approximately 1.6 kb.   |
| <b>Cloning Vector:</b>          | pCDNAI/amp (Invitrogen).  |
| <b>Description:</b>             | <p>Chemokine receptor expression vector, CCR1. The cDNA encoding the CCR1 was amplified from activated PBMC RNA using primers hybridizing to the 5' and 3' untranslated regions of CCR1. The PCR product was digested with BamHI and SalI and cloned into the BamHI and XhoI sites of pCDNAI/amp (XhoI and SalI are compatible ends for ligation). Constructs are selected for by adding Ampicillin (100 µg/mL) to the growth medium. Chloramphenicol (100 µg/ml) should be added 4.5-5.5 hours (OD<sub>600</sub> = 1.0) after seeding the culture. The plasmid can be transformed into E. coli strain HB101.</p> <p><u>CCR-1 primer sequences:</u><br/>CCR-1.5'BamHI (34mer): 5'-GCTCAGGATCCGCCAGAAACAAAGACTTCACGG-3'.<br/>CCR-1.3'SalI (34mer): 5'-CGATCGGTCGACGTTCTATGTTCCCCAGGATTCC-3'.</p> |
| <b>Special Characteristics:</b> | <p><a href="#">Contributor Provided Plasmid Map</a></p> <p><a href="#">Sequence file lot 97058</a></p> <p>This reagent is currently being provided as dried purified DNA stabilized in DNastable <i>PLUS</i>. Please see the notice for additional information and the protocol for reconstitution of dried DNA reagents. <a href="#">Dried DNA Notice</a></p>  |

---

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

**Recommended Storage:** Keep the reagent at room temperature in a dry storage cabinet or in a moisture barrier bag.

**Contributor:** Dr. Nathaniel Landau, Aaron Diamond AIDS Research Center, The Rockefeller University.

**References:** Morgenstern JP, Land H. Advanced mammalian gene transfer: high titre retroviral vectors with multiple drug selection markers and a complementary helper-free packaging cell line. *Nucleic Acids Res* **18**:3587-3596, 1990.

Deng H, Liu R, Ellmeier W, Choe S, Unutmaz D, Burkhart M, Di Marzio P, Marmon S, Sutton RE, Hill CM, Davis CB, Peiper SC, Schall TJ, Littman DR, Landau NR. Identification of a major co-receptor for primary isolates of HIV-1. *Nature* **381**:661-666, 1996.

**NOTE:** Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: Human CCR1 Expression Vector (pcCCR1) (Cat# 3321) from Dr. Nathaniel Landau." Also include the references cited above in any publications.

**Patent-pending. Requests from commercial organizations must be directed to the New York University Office of Industrial Liaison at the following email address: [sadhana.chitale@nyumc.org](mailto:sadhana.chitale@nyumc.org).**

**Last Updated:** November 09, 2018

---

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.