



NIH AIDS Reagent Program

20301 Century Boulevard
Building 6, Suite 200
Germantown, MD 20874
USA

Phone: 240 686 4740
Fax: 301 515 4015
aidsreagent.org

DATA SHEET

Reagent:	Human CCR2B Expression Vector (pcCCR2B)
Catalog Number:	3322
Lot Number:	150365
Release Category:	C
Provided:	5 µg of dried purified DNA stabilized in DNASTable <i>PLUS</i>
Cloning Site:	BamHI/XhoI cloning site The size of the insert is approximately 1083 bp.
Cloning Vector:	pcDNAI/amp (Invitrogen) Ampicilling resistant
Description:	Chemokine receptor expression vector, CCR2B.
Special Characteristics:	<p>The size of this construct is 5894 bp including the insert.</p> <p>The cDNA encoding CCR2B was amplified from activated PBMC RNA using primers hybridizing to the 5' and 3' untranslated regions of CCR2B. The PCR product was digested with BamHI and SalI and cloned into the BamHI and XhoI sites of pcDNAI/amp (SalI and XhoI are compatible ends for ligation).</p> <p>Constructs are selected for by adding Ampicillin (100 µg/mL) to the growth medium. Chloramphenicol (100 µg/ml) should be added 4.5–5.5 hours (OD₆₀₀ = 1.0) after seeding the culture.</p> <p><u>CCR-2b primer sequences:</u> CCR-2B.5'BamHI (35mer): 5'-GCTCAGGATCCTGAGACAAGCCACAAGCTGAACAG-3'. CCR-2B.3'XhoI (34mer): 5'-GTGCCTCTAGACTGAATGCGTGAGCCCTTGCTC-3'.</p> <p>Contributor provided plasmid map</p> <p>Plasmid map and sequence file lot 150365</p>

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

This reagent is currently being provided as dried purified DNA stabilized in DNastable *PLUS*. Please see the notice for additional information and the protocol for reconstitution of dried DNA reagents. [Dried DNA Notice](#)

- Recommended Storage:** Keep the reagent at room temperature in a dry storage cabinet or in a moisture barrier bag.
- Contributor:** Dr. Nathaniel Landau, Aaron Diamond AIDS Research Center, The Rockefeller University.
- References:** Morgenstern JP, Land H. Advanced mammalian gene transfer: high titre retroviral vectors with multiple drug selection markers and a complementary helper-free packaging cell line. *Nucleic Acids Res* **18**:3587-3596, 1990.
- Deng H, Liu R, Ellmeier W, Choe S, Unutmaz D, Burkhardt M, Di Marzio P, Marmon S, Sutton RE, Hill CM, Davis CB, Peiper SC, Schall TJ, Littman DR, Landau NR. Identification of a major co-receptor for primary isolates of HIV-1. *Nature* **381**:661-666, 1996.
- NOTE:** Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: Human CCR52B Expression Vector (pcCCR2B) (cat #3322) from Dr. Nathaniel Landau." Also include the references cited above in any publications.
- Patent-pending. Requests from commercial organizations must be directed to the New York University Office of Industrial Liaison at the following email address: sadhana.chitale@nyumc.org.**
- Last Updated:** September 19, 2017

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.