

NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

DATA SHEET

Reagent: Anti-HIV-1 gp120 Hybridoma (902)

Catalog Number: 521

Lot Number: 050754

С Release Category:

Characteristics:

Provided: 7×10^6 cells/vial. Viability is 60%.

Propagation RPMI 1640 supplemented with 200 U/ml penicillin, 2 mM L-glutamine, 0.05 mM Medium:

β-mercaptoethanol, 90%; fetal bovine serum, 10%.

Freeze Medium: Propagation medium, 40%; fetal bovine serum, 50%; DMSO, 10%.

Growth Important: Please see cell culture propagation instructions attached to data sheet. Cells **Characteristics:**

grow in suspension with some adherent cells. Growth is rapid. Doubling time is 24 hours.

Split every 3-4 days.

[BIO.A(2R) x A.BY]F1 mouse spleen cells x NS-1.

Description Thawing and Propagation of Hybridomas

Produces an IgG $_1$ κ monoclonal antibody (<u>Catalog #522</u>) which reacts with HIV-1_{LAV}/HTLV-III $_B$ gp120 in immunoprecipitation and Western blot assays. Neutralizes Special

infectivity and reacts with cells expressing the envelope proteins of these strains. Binds to the immunodominant hypervariable loop of gp120 (unpublished results). Does not recognize other HIV-1 strains. Mice were immunized via tail scratch with a vaccinia virus expressing gp160. Three-five weeks later, mice received two intraperitoneal inoculations

with the same virus. Three days after the second inoculation, mouse spleen cells were

fused with NS-1 cells.

Sterility: Negative for bacteria, fungi, and mycoplasma.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 09/20/2016 Page 1 of 2 Recommended Storage:

Liquid nitrogen.

Contributor: Dr. Bruce Chesebro.

References: Chesebro B, Wehrly K. Development of a sensitive quantitative focal assay for human

immunodeficiency virus infectivity. J Virol 62:3779-3788, 1988.

Pincus S, Wehrly K, Chesebro B. Treatment of HIV tissue culture infection with monoclonal

antibody-ricin A chain conjugates. J Immunol 142:3070-3075, 1989.

NOTE: Acknowledgment for publications should read "The following reagent was obtained

through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: Anti-HIV-1 gp120 Hybridoma (902) from Dr. Bruce Chesebro." Also include the references cited above in

any publications.

Scientists at for-profit institutions or who intend commercial use of this reagent

must contact the NIH Office of Technology Transfer,

Email: <u>NIAIDAIDSReagent@niaid.nih.gov</u>, before the reagent can be released. Please specify the name and a description of the intended use of the reagent.

Last Updated September 20, 2016

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 09/20/2016 Page 2 of 2