



**NIH AIDS Reagent Program**  
20301 Century Boulevard  
Building 6, Suite 200  
Germantown, MD 20874  
USA

Phone: 240 686 4740  
Fax: 301 515 4015  
aidsreagent.org

## DATA SHEET

**Reagent:** SIVmac239 Gag Peptide Set

**Catalog Number:** 6204

**Lot Number:** 009

**Provided:** 125 vials, 1.0 mg each, lyophilized.

**Sequence:** GenBank: M33262

**Purity:** Most peptides are >80% pure.

**Solubility:** See attached table. Peptides that are difficult to solubilize can almost always be dissolved in DMSO. Once a peptide is in solution, the DMSO can be slowly diluted with aqueous medium. Care must be taken to ensure that the peptide does not begin to precipitate out of solution.

**Please note: Cys and Met containing peptide sequences are highly unstable when stored in DMSO for prolonged periods of time; peptide decay is often observed within days to weeks. If DMSO is required, we recommend freshly-prepared working solutions from freeze dried aliquots of these peptides, and prepared pools.**

**Special Characteristics:** [Click here to see the peptides in this complete set.](#)

**Contributor:** DAIDS, NIAID.

**NOTE:** Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: SIVmac239 Gag Peptides - Complete Set."

**Last Updated:** March 09, 2015

---

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.