



## NIH AIDS Reagent Program

20301 Century Boulevard  
Building 6, Suite 200  
Germantown, MD 20874  
USA

Phone: 240 686 4740  
Fax: 301 515 4015  
aidsreagent.org

### DATA SHEET

**Reagent:** Anti-HIV-1 HXB2 IN 1-16 aa Polyclonal

**Catalog Number:** 756

**Lot Number:** 150042

**Provided:** 250 µl undiluted rabbit antiserum.

**Titer:** 1:20-50 (WB)

**Special Characteristics:** Antibodies were raised against HXB2 synthetic peptides. The antisera can be used to detect integrase in immunoprecipitation (IPT) and Western blot (WB) assays.  
  
Cat# 756 - aa 1-16: FLDGIDKAQDEHEKYH  
Cat# 757 - aa 23-34: ASDFNLPVAVAK  
Cat# 758 - aa 276-288: AGDDCVASRQDED

**Recommended Storage:** Keep at 4°C for short term storage and -80°C for long term storage. Avoid freeze-thaw cycles as reagent degradation may result.

**Contributor:** Dr. Duane P. Grandgenett.

**References:** Grandgenett DP, Goodarzi G. Folding of the multidomain human immunodeficiency virus type-I integrase. *Prot Sci* **3**:888-897, 1994.  
  
Bukrinsky MI, Sharova N, McDonald TL, Pushkarskaya T, Tarpley WG, Stevenson M. Association of integrase, matrix, and reverse transcriptase antigens of human immunodeficiency virus type 1 with viral nucleic acids following acute infection. *Proc Natl Acad Sci USA* **90**:6125-6129, 1993.

**NOTE:** Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: Anti-HIV-1 HXB2 IN 1-16 aa Polyclonal (cat# 756) from Dr. Duane P. Grandgenett." Also include the first reference listed above in any publications.

---

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.