

## NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

## DATA SHEET

HIV-1 CH293a Infectious Molecular Clone Reagent:

Catalog Number: 13539

Lot Number: 190378

С Release Category:

Provided: 5 μg of dried purified DNA stabilized in DNAstable PLUS

**Cloning Vector:** pUC57

Ampicillin resistant

**Cloning Site:** EcoRI/HindIII cloning site

The size of the insert is approximately 10,000 bp.

GenBank: KC156216

**Host Strain:** Plasmids can be propagated in STBL2 cells and grown at 37°C. Larger plasmids may

benefit from growth at 30°C. This construct may also be grown in other competent cells.

**Description:** A full length replication competent, infectious HIV-1 subtype C CH293a molecular clone.

The virus produced by this clone utilizes the co-receptor CCR5.

Special Characteristics:

This construct is approximately 14,000 bp including the insert.

Viral RNA was extracted from plasma samples from an HIV-1 infected female from Malawi. Consensus sequences were generated with clonally expanded viral sequences in the chronic infection samples. The sequences were cloned and combined within a single

plasmid vector (pBR322 or pUC57) as a complete proviral genome.

This reagent is currently being provided as dried purified DNA stabilized in DNAstable PLUS. Please see the notice for additional information and the protocol for reconstitution

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 12/22/2020 Page 1 of 2

## or aried DINA reagents. Dried DINA NOTICE

Recommended Storage:

Keep the reagent at room temperature in a dry storage cabinet or in a moisture barrier

bag.

Contributor: Dr. Feng Gao

References: Parrish, N. F., Gao, F., Li, H., Giorgi, E. E., Barbian, H. J., Parrish, E. H., Zajic, L., Iyer, S.

Parrish, N. F., Gao, F., Ll, H., Glorgi, E. E., Bardian, H. J., Parrish, E. H., Zajic, L., Lyer, S. S., Decker, J. M., Kumar, A., Hora, B., Berg, A., Cai, F., Hopper, J., Denny, T. N., Ding, H., Ochsenbauer, C., Kappes, J. C., Galimidi, R. P., West, A. P., Jr., Bjorkman, P. J., Wilen, C. B., Doms, R. W., O'Brien, M., Bhardwaj, N., Borrow, P., Haynes, B. F., Muldoon, M., Theiler, J. P., Korber, B., Shaw, G. M. and Hahn, B. H. (2013). Phenotypic properties of transmitted founder HIV-1. Proc Natl Acad Sci U S A, 110(17), 6626-33.

doi:10.1073/pnas.1304288110 PUBMED

NOTE: Acknowledgment for publications should read "The following reagent was obtained

through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: HIV-1 CH293a Infectious Molecular Clone from Dr. Feng Gao (cat# 13539)." Also include the reference

cited above in any publications.

Scientists at for-profit institutions or who intend commercial use of this reagent must contact the Duke University-Office of Research Contracts at the following email addresses: <u>david.kordys@duke.edu</u>or <u>bilyana.georgieva@duke.edu</u>, before

the reagent can be released.

**Last Updated:** December 22, 2020

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 12/22/2020 Page 2 of 2