

Streptococcus pneumoniae, Strain OREP3

Catalog No. NR-51850

Product Description:

Streptococcus pneumoniae (*S. pneumoniae*), strain OREP3 was derived from human wild-type *S. pneumoniae*, strain Wu2 (serotype 3) by natural selection using increasing concentrations of optochin. NR-51850 was produced by inoculation of BEI Resources seed lot 70029622 into Todd-Hewitt broth containing 0.5% (w/v) yeast extract and grown for 1 day at 37°C in an aerobic atmosphere with 5% CO₂. Broth inoculum was added to Tryptic Soy agar with 5% defibrinated sheep blood kolles, which were grown for 1 day at 37°C in an aerobic atmosphere with 5% CO₂ to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70062241

Manufacturing Date: 27JUL2023

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Hemolysis Motility (wet mount) Catalase VITEK® MS (MALDI-TOF)	Gram-positive cocci Report results α-hemolytic Report results Report results <i>S. pneumoniae</i>	Gram-positive cocci Circular, low convex, entire, mucoid and gray (Figure 1) α-hemolytic Non-motile Negative <i>S. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile¹ Thermo Scientific™ Oxoid Optochin Disc 1 day at 35°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood Optochin (Thermo Scientific™ DD0001T)	Report results	Resistant (0 mm)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1490 base pairs)	≥ 99% sequence identity to <i>S. pneumoniae</i> type strain (GenBank: NR_028665.1)	99.9% sequence identity to <i>S. pneumoniae</i> type strain (GenBank: NR_028665.1)
Purity (post-freeze) 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability (post-freeze)	Growth	Growth

¹Zone of Inhibition; Interpretation Guideline: Burton, R. L. and M. H. Nahm. "Development of a Fourfold Multiplexed Opsonophagocytosis Assay for Pneumococcal Antibodies against Additional Serotypes and Discovery of Serological Subtypes in *Streptococcus pneumoniae* Serotype 20." *Clin. Vaccine Immunol.* 19 (2012): 835-841. PubMed: 22518015.

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

20 OCT 2023

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

