

Certificate of Analysis for NR-44210

Complementary DNA from Biomphalaria glabrata, Strain NMRI

Catalog No. NR-44210

This reagent is the tangible property of the U.S. Government.

Product Description: Complementary DNA (cDNA) was synthesized from RNA extracted from *Biomphalaria glabrata* (*B. glabrata*), strain NMRI.

Lot¹⁻³: 61895539 Manufacturing Date: 21MAY2013


TEST	SPECIFICATIONS	RESULTS
Concentration	Report results	0.5 μg in 20 μL per vial (0.025 μg/μL)
OD ₂₆₀ /OD ₂₈₀ Ratio	1.80 to 2.00	1.98
Qualification by RT-PCR Amplification of myoglobin gene ⁴	~ 190 base pair amplicon	~ 190 base pair amplicon (Figure 1)

QC testing was performed by the Biomedical Research Institute, Rockville, MD (NIH-NIAID Contract HHSN272201000005I)

Figure 1: Amplification of Myoglobin Gene for Qualification by RT-PCR

Lane 1: 100 base pair ladder

Lane 2: 190 base pair amplicon from *B. glabrata* myoglobin gene



Date: 18 OCT 2013

Signature:

Title: Technical Manager, BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by the contributor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources

www.beiresources.org

Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898

²Snails were treated with 100 µg/mL ampicillin overnight to reduce rotifer and bacterial contamination

³cDNA was amplified from total RNA by SMART™ cDNA Synthesis Kit (Clontech) according to the manufacturer's instructions.

⁴Primers were designed to amplify the nucleotide region 108 to 296 of B. glabrata myoglobin gene (GenBank AY314976).