

***Peptoclostridium difficile*, Strain CD131**

**Catalog No. NR-43617**

**Product Description:** *Peptoclostridium difficile* (*P. difficile*; also referred to as *Clostridium difficile*), strain CD131 is a toxigenic strain isolated in March 2010 from the stool of a human patient diagnosed with a recurrent *Clostridium difficile* infection in Ann Arbor, Michigan, USA.

**Lot<sup>1</sup>: 63312215**

**Manufacturing Date: 13FEB2015**

TEST	SPECIFICATIONS	RESULTS
<b>Phenotypic Analysis</b> Cellular morphology Colony morphology <sup>2</sup>  Hemolysis <sup>2</sup> Motility (wet mount) Biochemical tests: Esculin hydrolysis Gelatin hydrolysis VITEK <sup>®</sup> MS (MALDI-TOF)	Gram-positive rods Report results  Report results Report results  Positive Positive Consistent with <i>P. difficile</i>	Gram-positive rods Irregular, flat, undulate and gray (Figure 1) Non-hemolytic Motile  Positive Positive Consistent with <i>P. difficile</i>
<b>Genotypic Analysis</b> Sequencing of 16S ribosomal RNA gene (~ 1450 base pairs)	Consistent with <i>P. difficile</i>	Consistent with <i>P. difficile</i> <sup>3</sup>
<b>PCR Assay of Extracted DNA</b> Presence of <i>P. difficile</i> -specific genes <sup>4</sup> Triose phosphate isomerase ( <i>tpi</i> ) Presence of toxin genes <sup>4</sup> Toxin A ( <i>tcdA</i> ) Toxin B ( <i>tcdB</i> )	~ 230 base pairs amplicon  ~ 369 base pairs amplicon ~ 160 base pairs amplicon	~ 230 base pairs amplicon  ~ 369 base pairs amplicon ~ 160 base pairs amplicon
<b>Purity (post-freeze)</b> Anaerobic growth <sup>5</sup> Aerobic growth <sup>6</sup>	Growth consistent with <i>P. difficile</i> No growth	Growth consistent with <i>P. difficile</i> No growth
<b>Viability (post-freeze)<sup>2</sup></b>	Growth	Growth

<sup>1</sup>NR-43617 was produced by inoculation of the deposited material into Modified Reinforced Clostridial medium and incubated for 24 hours at 37°C in an anaerobic atmosphere (< 0.5% O<sub>2</sub>; Remel™ Anaero Pack-Anaero™ R681001). The material from the initial growth was passaged once in Modified Reinforced Clostridial medium for 24 hours under propagation conditions to produce this lot.

<sup>2</sup>23 hours on Tryptic Soy agar with 5% defibrinated sheep blood under propagation conditions


<sup>3</sup>≥ 99.4% identical to GenBank: AVHQ01000046.1 (*P. difficile*, strain CD131)

<sup>4</sup>Lemee, L., et al. "Multiplex PCR Targeting *tpi* (Triose Phosphate Isomerase), *tcdA* (Toxin A), and *tcdB* (Toxin B) Genes for Toxigenic Culture of *Clostridium difficile*." *J. Clin. Microbiol.* 42 (2004): 5710-5714. PubMed: 15583303.

<sup>5</sup>Purity of this lot was assessed for 7 days on Tryptic Soy agar with 5% defibrinated sheep blood under propagation conditions.

<sup>6</sup>Purity of this lot was assessed for 7 days on Tryptic Soy agar with 5% defibrinated sheep blood in an aerobic atmosphere with 5% CO<sub>2</sub>.

**Figure 1**



Date: 21 MAY 2015

Signature: 

BEI Resources Authentication

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection.  
You are authorized to use this product for research use only. It is not intended for human use.

