

***Pseudomonas aeruginosa*, Strain MRSN 16345**

**Catalog No. NR-51582**

This reagent is the tangible property of the U.S. Government.

**Product Description:**

*Pseudomonas aeruginosa* (*P. aeruginosa*), strain MRSN 16345 was isolated in 2013 from human urine as part of a surveillance program in the United States. *P. aeruginosa*, strain MRSN 16345 was deposited as sensitive to amikacin, gentamicin, meropenem and tobramycin, intermediately resistant to cefepime and imipenem and resistant to aztreonam, ceftazidime, ciprofloxacin, levofloxacin and piperacillin/tazobactam.

**Lot: 70025070<sup>1</sup>**

**Manufacturing Date: 27JUN2019**

TEST	SPECIFICATIONS	RESULTS
<b>Phenotypic Analysis</b> Cellular morphology Colony morphology <sup>2</sup>  Motility (wet mount) VITEK <sup>®</sup> 2 (GN card)	Gram-negative rods Report results  Report results <i>P. aeruginosa</i> (≥ 89%)	Gram-negative rods Irregular, flat, undulate, smooth and green (Figure 1) Motile <i>P. aeruginosa</i> (98%)
<b>Antibiotic Susceptibility Profile<sup>3</sup></b> VITEK <sup>®</sup> (AST-GN81 Card) Ampicillin Amoxicillin/clavulanic acid Piperacillin/tazobactam Cefazolin Cefoxitin Ceftazidime Ceftriaxone Cefepime Meropenem Amikacin Gentamicin Tobramycin Ciprofloxacin Levofloxacin Tetracycline Nitrofurantoin Trimethoprim/sulfamethoxazole	Report results Report results Resistant Report results Report results Resistant Report results Intermediate Sensitive Sensitive Sensitive Sensitive Resistant Resistant Report results Report results Report results	Resistant (≥ 32 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (32 µg/mL) <sup>4</sup> Sensitive (1 µg/mL) Sensitive (≤ 2 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Resistant (≥ 4 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (256 µg/mL) ≥ 320 µg/mL <sup>5</sup>
<b>Genotypic Analysis</b> Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>P. aeruginosa</i> , strain MRSN 16345 (GenBank: RXVR01000117.1)	100% sequence identity to <i>P. aeruginosa</i> , strain MRSN 16345 (GenBank: RXVR01000117.1)
<b>Purity (post-freeze)<sup>6</sup></b>	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
<b>Viability (post-freeze)<sup>2</sup></b>	Growth	Growth

<sup>1</sup>NR-51582 was produced by inoculation of the depositor material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot.

<sup>2</sup>1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar


<sup>3</sup>Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

<sup>4</sup>Susceptibility results for this antibiotic are within one doubling dilution of specification, which is considered an equivalent result.

<sup>5</sup>Trimethoprim/sulfamethoxazole MIC interpretive standards are not available for *P. aeruginosa*, however most clinical isolates are resistant to trimethoprim/sulfamethoxazole. For more information, please refer to Köhler, T., et al. "Multidrug Efflux in Intrinsic Resistance to Trimethoprim and Sulfamethoxazole in *Pseudomonas aeruginosa*." *Antimicrob. Agents Chemother.* 40 (1996): 2288-2290. PubMed: 9036831.

<sup>6</sup>Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere with and without 5% CO<sub>2</sub> on Tryptic Soy agar.

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch

11 DEC 2019

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

