

Certificate of Analysis for NR-59584

Genomic RNA from Influenza A Virus, A/Baltimore/JH-0440/2022 (H3N2)

Catalog No. NR-59584

Product Description:

Genomic RNA was isolated from a preparation of cell lysate and supernatant from Madin-Darby canine kidney SIAT-1 (MDCK-SIAT1) cells infected with influenza A virus, A/Baltimore/JH-0440/2022 (H3N2) using QIAamp® Viral RNA Mini Kit (Qiagen® 52906). The viral genomic RNA is in a background of cellular nucleic acid and carrier RNA.

Lot: 70063751 Manufacturing Date: 18DEC2023

TEST	SPECIFICATIONS	RESULTS
Genotypic Analysis Sequencing of species-specific region (~ 870 nucleotides)	Consistent with source virus	Consistent with source virus ¹
Functional Activity by RT-PCR Amplification ² Matrix gene	~ 1000 base pair amplicon	~ 1000 base pair amplicon
Estimated Concentration (post-dilution) by Qubit [®] Measurement (Viral, Cellular and Carrier) ³	Report results	5.3 ng per 100 μL (0.053 μg/mL)
Estimated Amount per Vial ³	Report results	5.3 ng
Genome Copy Number Using BioRad QX200 Droplet Digital PCR (ddPCR™) System (Post vial; 12 replicates)	Report results	2.6 × 10 ⁸ NDU/mL ⁴
Virus Inactivation 10% of total yield inoculated on MDCK-SIAT1 cells and evaluated for cytopathic effect and HA after serial passage ⁵	No viable virus detected	No viable virus detected

Sequence information for influenza virus A/Baltimore/JH-0440/2022 (H3N2) is not available in the NCBI database; nucleotide sequence obtained for NR-59584 lot 70063751 is identical to the source virus.

/Sonia Bjorum Brower/ Sonia Bjorum Brower

<u>09 APR 2024</u>

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898

²Amplified using iTaq[™] Universal SYBR Green One-step Kit (Bio-Rad[®] 172-5151) with 5 μL of NR-59584 in a 50 μL reaction

³Measurement is determined pre-vial prior to dilution due to the limit of detection of the quantification method

⁴NDU; NAAT-detectable units

⁵Use of the QlAamp® Viral RNA Mini Kit has been demonstrated to consistently inactivate 100% of influenza A viruses as shown by the absence of cytopathic effect (CPE) and HA after plating the entire extract on virus-susceptible cells for two passages.